



# DSX 系列

## 手提式 压力蒸汽灭菌器 (医疗) 手提式 高压蒸汽灭菌器 (非医疗)

Portable Pressure Steam Sterilizer



型号	DSX-280KB24 (医疗) DSX-24L (非医疗)	DSX-280KB30 (医疗) DSX-30L (非医疗)
容积	24 立升	30 立升

## Product Features 产品特点

### 实行微电脑智能控制系统：

采用“LED”数显运行工作循环程序，灭菌结束（报警）后自动停机；

### 安全联锁装置：

采用电子与机械互动的安全 构，确保锅内有压力时联锁装置自动锁紧锅盖，避免了误操作或使用不当而产生的不安全因素；

### 锅盖启闭装置：

采用移位式快开盖结构，使锅盖开启与锁紧安全可靠、灵活轻巧；

### 预置程序：

可预置固定程序对灭菌物品的灭菌；

### 记忆存储系统：

可存储上次灭菌程序，便于二次操作；

### 操控台动态：

数显屏与指示灯可显示设定的灭菌工作动态与运行的功能提示；

### 安全配置：

锁盖联动装置、电控式超温过压保护系统、防干烧保护功能、机械式安全泄压阀、漏电保护装置。

### Intelligent microcomputer control system:

Using LED digital control system; After sterilization finished, sterilizer will automatic alarm and shut down.

### Interlock device:

The control system can real time supervision the sterilizer chamber's temperature and pressure. Only under the safety temperature and pressure can open the sterilizer's cover. Make sure the operator's safety.

### The structure of cover's opening:

Using shift place type of quickopen cover to improve sealing and safety performance.

### Preset program:

Can preset fixed program that sterilize solid body.

### Memory storage system:

Store the last operation program so convenient next operation.

### Digital display screen:

Digital display screen and working indicator light display running state.

### Safety devices:

Interlock device, overtemperature and over pressure protecting system, burned protection, safety relief valve, current leakage protector and automatic fault detection.

## Technical Parameter 技术参数

参数 / 型号 Parameter / Model	DSX-280KB24 (医疗) DSX-24L (非医疗)	DSX-280KB30 (医疗) DSX-30L (非医疗)
配置 Configuration	自动控制 Automatic Control	自动控制 Automatic Control
材质 / 锅体壁厚 (mm) Material / Pot wall thickness(mm)	全不锈钢 / 2.0 Stainless steel / 2.0	全不锈钢 / 2.0 Stainless steel / 2.0
容积 / 毛重 Volume / GW	24L / 25Kg	30L / 27.5Kg
电压 / 功率 Voltage / Power	220V / 2Kw	220V / 2Kw
包装尺寸 (长 x 宽 x 高 mm) Packing size (L X W X H mm)	500 X 500 X 800	500 X 500 X 900
灭菌室容积 (直径 x 高度 mm) Sterilization chamber volume (D X H mm)	280 X 495	280 X 585
时间设定范围 Time setting range	0-999min	0-999min
额定工作压力 Rated working pressure	0.142MPa	0.142MPa
灭菌温度设定范围 Sterilization temperature setting range	50~126℃	50~126℃